

InVite™

liquid lure

**ATTRACTS YELLOW JACKETS, EUROPEAN
HORNETS AND FRUIT FLIES**

Concentrated Formula • Contains No Pesticides

Directions For Use

Shake well and pour 1-4 oz. (depending upon level of alternative food in the area) of InVite Liquid Lure into the cup of the PFT Station and add water, soda or beer to fill the cup about 1/2 full. This product will also attract phorid flies and can be enhanced for this use by adding a small spoonful of Rockwell's InVite Fly Lure to the mixture. Hang the PFT Station in areas where target pests are present. To create a killing station, spray or wipe the PFT with an odorless residual insecticide. Replace when the cup is full or after one month. Submerge Station in water prior to disposing of cup to kill any insects left alive. For use in ILTs, place InVite Liquid Lure in a small container in or on the ILT. This product may also be used in other insect traps. To use on a glueboard, place a larger cottonball on the glue and saturate with InVite Liquid Lure. For greater convenience on glueboards, use InVite Fruit Fly Lure (gel) instead. Avoid ingestion and inhalation of fumes for extended periods. Store out of reach of children. For long term storage, avoid excessive heat. Used or leftover InVite Liquid Lure and dead insects may be disposed of on site.

Net Contents: 16 fl. oz.

 **Rockwell Labs Ltd**
creating the future of pest control

www.rockwelllabs.com
toll free: 866 788 4101
EPA Est. No. 73079-MO-1